

RadiForce GS520

Eizo RadiForce GS520

21.3" Monochrome LCD Monitor

Featuring high-resolution and high-definition, the RadiForce GS520 is designed specifically for displaying digital mammography images with a choice of distinctive glare panel to meet diversified environment usage.



High-Definition and High-Density with Minute Pixels

0.165mm pixel size is the smallest pixel size among current EIZO medical monitors providing high-definition and highdensity display without shadings appearing granular.

12-Bit Simultaneous Grayscale Display

Along with EIZO's new frame rate control technology (patent pending), 4,096 (12-bit) grayscale tones can be displayed simultaneously from an abundant palette of 13,771 tones for high-definition digital mammography.

Distinctive Glare Panel

Many current monitors use anti-glare panels, which have a waffled surface to diffuse the reflection of surrounding light. However, this also diffuses the light from the monitor's backlight, which affects the outline of an image, causing a viewing strain on the observer. With a distinctive glare panel, there is no diffused reflection by the waffled surface, and the outline of the image is displayed more clearly. Therefore, the subtle images such as calcification and masses will be easier to observe reducing the strain on the observer.

Sharper Resolution for Clearer Image Outline

MTF (Modulation Transfer Function) is a parameter to define the display sharpness and resolution of a monitor panel. The glare panel has better MTF characteristics

than the anti-glare panel meaning the glare panel can be expected to display the outline of a detailed image more sharply.

DUE for Brightness Uniformity

The Digital Uniformity Equalizer (DUE) function provides optimum backlight luminance uniformity which is considered difficult to attain due to the characteristics of LCD monitors.

Calibration Mode Selection

Selectable with the front panel buttons, the CAL Switch function allows for various calibration modes of different modalities such as digital mammography, ultrasound, and MRI images. Furthermore, with ScreenManager Pro for Medical installed, auto mode settings can be made with the Auto CAL Switch function.

Brightness Stabilization

At startup or upon wakeup, the EIZO patented drift correction function quickly stabilizes the brightness level. In addition, a sensor measures the backlight brightness and compensates for brightness fluctuations caused by the ambient temperature and the passage of time.

Wide Viewing Angles

Viewing angles are 170° wide with minimal color shift when viewed from the side.



RadiForce GS520

Backlight Saver

With ScreenManager Pro for Medical utility software installed, the Backlight Saver function allows for the monitor's backlight to turn off when the screen saver is activated and the monitor's backlight to turn on when the computer comes out of the screen saver mode. This function helps to reduce power consumption when the monitor is used for a prolonged period of time.

Ergonomic Design

Narrow bezels require less eye movement in a multipanel environment. The height and pivot adjustable stand offers added flexibility and viewing comfort.

Compatibility with Various Graphics Boards

Supports various input signals for compatibility with a wide range of graphics boards so the user can select the board according to budget and performance requirements. This includes single link and dual link signals, as well as packed-pixel signal for users requiring 12-bit or 10-bit simultaneous display.

Pairing Monitors

EIZO offers the same color and grayscale for two monitors bundled as a pair in one box.

Clear Base and Blue Base Backlights

Available with either a clear base or blue base backlight.

Full 5-Year Warranty

EIZO and its authorized distributors offer a five-year limited warranty.

Brightness Stability Within Usage Time Guaranteed

EIZO's confidence in its product quality extends to brightness stability which is also covered during the usage time specified in the warranty.

Customer Assurance

Meets the strictest medical, safety and EMC emissions standards including TÜV/GM, UL60601-1, CE, CB, CSA C22.2 No.601-1, VCCI, FCC, and FDA 510(k).

ISO 13485 Certification

EIZO has acquired ISO 13485 certification for its quality management system. Approval to this standard demonstrates EIZO's ability to consistently meet customer requirements for these devices and services.





RadiForce GS520

Specifications

Model Variations GS520-CLG: Clear Base with Glare Panel

GS520-BLG: Blue Base with Glare Panel

GS520-CL: Clear Base GS520-BL: Blue Base

Cabinet Color Black

Panel TFT Monochrome LCD Panel (IPS)

Active Display Size (H x V)
Viewable Image Size
Native Resolution

Pixel Pitch

337.9 x 422.4 mm
Diagonal: 540 mm
2048 x 2560
0.165 x 0.165 mm

Grayscale Tones 4,096 from a palette of 13,771 tones

500 cd/m2

Viewing Angles (H, V) 170°, 170°

Brightness 700 cd/m2 (typical)

Recommended Brightness

for Calibration

Contrast Ratio 800:1 (typical)
Response Time 100 ms (typical)

Scanning Frequency (H, V) 31-135 kHz, 19-51 Hz

Frame synchronous mode: 49 - 51 Hz

Dot Clock 300 MHz

Input Signals DVI Standard 1.0 Input Terminal DVI-D 24 pin

USB Ports 1 upstream, 2 downstream

USB Standard Rev. 2.0

Power Requirements AC 100 -120 V, 200-240 V: 50 / 60 Hz

Power Consumption 90 watts

Power Save Mode
Sensor
Power Management

Less than 1.2 watts
Backlight Sensor
DVI DMPM

OSD Languages English, French, German, Italian, Japanese, Simplified Chinese, Spanish, Swedish,

Traditional Chinese

Net Weight With Stand: 9.6 kg

Without Stand: 6.6 kg

Hole Spacing VESA standard (100 x 100 mm)

Certifications and Standards CE (Medical Device Directive), TÜV/GM (EN60601-1), cTÜVus (UL60601-1, CSA

C22.2 No. 601-1), CB (IEC60601-1), VCCI-B, FCC-B, Canadian ICES-003-B, c-Tick, FDA

510(k) for Mammography and General Radiography, RoHS

Supplied Accessories AC power cord, user's manual, signal cable (DVI-D ~ DVI-D), USB cable, EIZO LCD

Utility Disk (RadiCS LE, ScreenManager Pro for Medical software), ScreenCleaner

Tel: +31 (0) 182 379 377

Fax: +31 (0) 182 379 374

(GS520-CLG, GS520-BLG only), warranty card